LuxaCore Z Dual

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English

Product description

LuxaCore Z Dual is an automatic mixing, dual-curing composite that has been specially developed for all types of core build-ups and build-up fillings. The curing time can be determined by the user due to the light-curing option. It is possible to apply LuxaCore Z Dual directly using the Intraoral Tip and the Endo Tip. LuxaCore Z Dual is radiopaque.

Intended purpose

- Fabrication of all types of core build-ups
- Cementation of root posts

Limitations of use

Do not use directly on the exposed pulp.

Indication

Extensive loss of tooth substance due to caries, trauma and development disorders.

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist

Notes for use

- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.
- Light-curing units should emit at 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm².
- Place the light-curing unit as close as possible to the material.

Using the Automix Cartridge (for image, see »Handling«)

Use of the Automix Cartridge

1. Push up the lever [A] at the back of the dispensing gun and pull back the slide [B] completely.
2. Lift up the plastic lever [C], insert the Automix Cartridge and lock it using the plastic lever [C].

Note: Make sure that the notches on the Automix Cartridge line up with those on the dispensing unit.

Attaching a mixing tip

1. Remove the cartridge cap [D] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

Note: Make sure that the notches on the mixing tip line up with those on the Automix Cartridge.

3. Twist the mixing tip clockwise by 90° to lock it.
4. Attach an Intraoral Tip [F], if applicable.

Applying the material

- When using an Automix Cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- Mix the material in the mixing tip using the lever [G] on the dispensing gun and apply directly.

Note: After the application, leave the used mixing tip on the Automix Cartridge as a seal!

Removal of the Automix Cartridge

1. Push up the lever [A] at the back of the dispensing gun and pull back the slide [B] completely.
2. Lift up the plastic lever [C] and remove the cartridge.

Using the Smartmix Syringe (for image, see »Handling«)

1. Before attaching the mixing tip, remove the cap [A] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [B].

Note: In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Smartmix Syringe [C].

3. Twist the mixing tip clockwise by 90° to lock it.
4. Attach an Intraoral Tip or Endo Tip.
5. The material is mixed in the mixing tip when ejected and can then be applied directly.

Note: When using a Smartmix Syringe for the first time, squeeze out a pea-sized amount of material and discard it! After application, leave the used mixing tip on the Smartmix Syringe as a seal!

Disinfection of the Automix Dispenser

- The Automix Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear.
- Disinfection can be carried out using a commercially available immersion disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch-Institute). The Automix Dispenser is also autoclavable.

Recommended use

Carry out core build-ups

1. Isolate the tooth prior to preparation using a rubber dam.
2. Remove existing fillings and caries.
3. If the bonding agent used requires an additional etching step, then etch following the manufacturer's instructions for the etching gel used.

Note: The use of the LuxaBond Universal or the LuxaBond Total Etch dual-curing adhesive system is recommended. When using LuxaBond Total Etch, be sure to etch the affected enamel and dentin surfaces. Observe the manufacturer's instructions.

Caution! Unprotected pulp.

- Prior to bonding in deep cavities in areas close to pulp, provide protection, for example by using a small quantity of calcium hydroxide and a thin layer of glass ionomer cement.

4. Apply a bonding agent in accordance with the manufacturer's instructions.
5. In order to simplify the application of LuxaCore Z Dual, place a matrix around the prepared tooth and apply LuxaCore Z Dual directly into the cavity. If a pre-fabricated core build-up form is used, apply LuxaCore Z Dual directly into the core build-up form and place on the prepared tooth.

Note: LuxaCore Z Dual can be formed within the working time of 1:30 minutes using a standard instrument, e.g. a Heidemann spatula. In order to enhance processing, coat the instrument with a bonding agent as a separating medium.

6. Allow LuxaCore Z Dual to set intra-orally for approx. 5 minutes or light-cure with a suitable curing light for a minimum of 40 seconds (maximum layer thickness ≤ 4 mm). In the case of thinner layers (≤ 2 mm) exposure to the curing light for 20 seconds is sufficient.

Note: With chemical curing, lower temperatures (e.g. when working in the area of the front teeth) lead to longer setting times.

7. If a matrix is used, remove it only after the material has cured completely.
8. Prepare the desired restoration form using standard preparation instruments.
9. If a temporary resin-based restoration is to be applied to the prepared tooth, isolate the preparation using a suitable separation medium (e.g. Vaseline).

Cementation of root posts

1. Prepare the root post as per manufacturer's instructions.
2. Prepare the root canal in accordance with the requirements of the respective root post.
3. If the bonding agent used requires an additional etching step, then etch inside the canal following the manufacturer's instructions for the etching gel used.

Note: The use of the LuxaBond Universal or the LuxaBond Total Etch dual-curing adhesive system is recommended. When using LuxaBond Total Etch, be sure to etch the affected enamel and dentin surfaces. Observe the manufacturer's instructions.

4. Apply a bonding agent in accordance with the manufacturer's instructions. If an adhesive system with a light-curing component is used, use a paper tip to remove any excess from the canal prior to the exposure to light. Otherwise the lumen may be blocked.
5. Apply LuxaCore Z Dual in the canal lumen. When doing so keep the end of the tip constantly immersed in the material. If desired, it is additionally possible to coat the root post with LuxaCore Z Dual. The LuxaCore Z Dual working time is 1:30 minutes.

Note: When using LuxaBond Total Etch, the working time of LuxaCore Z Dual in the root canal is reduced to approx. 30 seconds due to the contact with the Pre-Bond. It is also possible that the working time may be reduced with other bonding agents. Observe the manufacturer's instructions.

6. Place the root post.
7. Light-cure the material with a suitable lamp (halogen or LED lamp) having a light intensity of ≥ 600 mW/cm² for 20 seconds.
8. Allow the material to self-cure for approx. 5 minutes. This guarantees complete curing, also in the areas that cannot be reached by light.

Note: With LuxaCore Z Dual it is possible to finish the coronal portion of the tooth as a core build-up. Due to the flowable consistency of the material, it is recommended that work be carried out using a matrix (see »Carry out core build-up«).

Timing

Working time (including mixing time)*	≈ 1:30 min
Chemical curing	Intraoral 5:00 min
Light-curing	20 seconds (layer ≤ 2 mm) 40 seconds (layer ≤ 4 mm)

* The specified times are valid at a room temperature of 23 °C and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times.

Interactions

Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Complies with ISO 4049:2019, type 2, class 3

Storage/disposal

- Store in a dry place at room temperature (2 - 25 °C/36 - 77 °F)!
- After application, leave the used mixing tip on the Automix Cartridge or Smartmix Syringe as a seal!
- Do not use after the expiration date.
- Disposal must comply with national regulations.

Composition

Dental glass, Silica (sil.), bis-GMA,UDMA, HDDMA, DDDMA, MMHE, additives.

Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

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